



BD FACS Canto II: 3-laser, 8-color (4-2-2)

Lasers		Channel Name	Filter		Other fluorochromes in this filter
Violet - 405nm	V1	Pacific Blue	450/40		BV421, AF405, Dapi
	V2	AmCyan	510/50	502 LP	BV510, Zombie Aqua
Blue - 488nm	FSC	SSC	488/10		
	B1	FITC	530/30	502 LP	GFP, AF488, CFSE
	B2	PE	585/42	556LP	tdTomato, PI, RFP
				610LP	
	B3	PerCP-Cy5.5	670LP	655 LP	PerCP
	B4	PE-Cy7	780/60	735 LP	
Red - 633nm	R1	APC	660/20		AF647
				685 LP	
	R2	APC-Cy7	780/60	735 LP	